

**CIL**Cambridge Isotope Laboratories, Inc.  
[www.isotope.com](http://www.isotope.com)

SOLVENTS

# LC/NMR Solvents

Available from CIL



## LC/NMR Solvents

Catalog No.	Compound	Size
DLM-9	Acetone-D <sub>6</sub> (D, 99.9%)	1 L
DLM-21	Acetonitrile-D <sub>3</sub> (D, 99.8%)	1 L
DLM-22	Acetonitrile-D <sub>3</sub> (D, 96-97%)	1 L
DLM-20	Acetonitrile-D <sub>3</sub> HPLC grade (D, 96-97%) + 15-20% D <sub>2</sub> O	1 L
DLM-4	Deuterium oxide (D, 99.9%)	1 kg
DLM-4-99	Deuterium oxide (D, 99%)	1 kg
DLM-6	Deuterium oxide "100%" (D, 99.96%)	1 kg
DLM-24	Methanol-D <sub>4</sub> (D, 99.8%)	1 L

All LC/NMR solvents can be packaged in convenient heavy-duty Pyrex bottles with GL45 thread to fit directly on the LC system.

## LC/NMR Reference Standards

Catalog No.	Compound	Size
DLM-76	1% Chloroform in acetone-D <sub>6</sub>	50 mL
DLM-5032	50 mM Sodium acetate plus 1% water in D <sub>2</sub> O	50 mL
DLM-7005	1% Water in D <sub>2</sub> O	50 mL

Please contact CIL to request a quote for pricing and availability.

To place an order please contact CIL: | t: 978.749.8000 | 1.800.322.1174 (North America) | [cilsales@isotope.com](mailto:cilsales@isotope.com)

For international inquiries, please contact our International Customer Service Department at [intlsales@isotope.com](mailto:intlsales@isotope.com).

NMR\_RSCH\_LCNMR  
2/14 Supersedes all previously published literature