

APO-14 Kit Specification Sheet for SIS Peptide Mixture (Human)

Protein	UniProt Acc. No.	Protein MW (Da)	Peptide	LLOQ (ng/mL)	ULOQ (ng/mL)
Apolipoprotein A-I	P02647	28078.62	ATEHLSTLSEK	36970	2957588
Apolipoprotein A-I	P02647	28078.62	VQPYLDDFQK	135073	2701451
Apolipoprotein A-II	P02652	8707.91	EQLTPLIK	9143	731464
Apolipoprotein A-II	P02652	8707.91	SPELQAEAK	5084	406694
Apolipoprotein A-IV	P06727	43402.53	LGEVNTYAGDLQK	27203	544067
Apolipoprotein A-IV	P06727	43402.53	LLPHANEVSQK	2848	227863
Apolipoprotein A-IV	P06727	43402.53	SLAPYAQDTQEK	3967	317332
Apolipoprotein B-100	P04114	512858.3	FPEVDLTK	83216	832160
Apolipoprotein C-I	P02654	6630.58	EFGNTLEDK	2153	43056
Apolipoprotein C-I	P02654	6630.58	TPDVSSALDK	1185	23705
Apolipoprotein C-II	P02655	8204.17	ESLSSYWESAK	4283	85652
Apolipoprotein C-II	P02655	8204.17	TAAQNLYEK	615	49225
Apolipoprotein C-IV	P55056	11815.57	AWFLESK	538	10755
Apolipoprotein D	P05090	19303.08	NILTSNNIDVK	19567	391339
Apolipoprotein D	P05090	19303.08	VLNQELR	2972	237753
Apolipoprotein E	P02649	34236.68	AQAWGER	2128	170203
Apolipoprotein E	P02649	34236.68	LGPLVEQGR	6044	120886
Apolipoprotein E	P02649	34236.68	SELEEQLTPVAEETR	10271	205420
Apolipoprotein E	P02649	34236.68	WELALGR	20542	410840
Apolipoprotein L1	O14791	41128.68	ALADGVQK	944	75512
Apolipoprotein L1	O14791	41128.68	ALDNLAR	267	21353
Apolipoprotein L1	O14791	41128.68	ANLQSVPHASASR	196	15658
Apolipoprotein M	O95445	21253.29	AFLTPR	2421	193693
Apolipoprotein(a)	P08519	499185.5	GTYSTTVTGR	8630	690388
Beta-2-glycoprotein 1	P02749	36254.6	ATVVYQGER	8605	688431
Clusterin	P10909	50062.56	IDSLENDR	6571	525657
Clusterin	P10909	50062.56	TLLSNLEAK	9093	181859

Cambridge Isotope Laboratories, Inc. is the exclusive distributor of MRM Proteomics Inc.

To place an order please contact CIL: | t: 978.749.8000 | 1.800.322.1174 (North America) | cilsales@isotope.com

For international inquiries, please contact our International Customer Service Department at intl@isotope.com.